



PAQR7 Antibody, Biotin conjugated

Product Code	CSB-PA768224LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WK9
Immunogen	Recombinant Human Membrane progesterin receptor alpha protein (1-72AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane progesterin receptor alpha (mPR alpha) (Membrane progesterone P4 receptor alpha) (Membrane progesterone receptor alpha) (Progesterone and adipoQ receptor family member 7) (Progesterin and adipoQ receptor family member 7) (Progesterin and adipoQ receptor family member VII), PAQR7, MRPA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PAQR7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.