



PLCZ1 Antibody, Biotin conjugated

Product Code	CSB-PA768240LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YW0
Immunogen	Recombinant Human 1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase zeta-1 protein (1-415AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase zeta-1 (EC 3.1.4.11) (Phosphoinositide phospholipase C-zeta-1) (Phospholipase C-zeta-1) (PLC-zeta-1) (Testis-development protein NYD-SP27), PLCZ1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	PLCZ1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.