





FAM122B Antibody

Product Code	CSB-PA768758LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z309
Immunogen	Recombinant Human Protein FAM122B protein (2-247AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM122B, FAM122B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM122B

Image



Western blot
All lanes: FAM122B antibody at 5µg/ml + Mouse brain tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 27, 28, 30, 20 kDa
Observed band size: 27 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA768758LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.