




LYSMD3 Antibody

Product Code	CSB-PA768760LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z3D4
Immunogen	Recombinant Human LysM and putative peptidoglycan-binding domain-containing protein 3 protein (1-127AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LysM and putative peptidoglycan-binding domain-containing protein 3, LYSMD3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LYSMD3

Image



Western blot
All lanes: LYSMD3 antibody at 6µg/ml + Mouse heart tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 35, 15, 18 kDa
Observed band size: 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.