




CYP2U1 Antibody

Product Code	CSB-PA768766LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z449
Immunogen	Recombinant Human Cytochrome P450 2U1 protein (363-496AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 2U1 (EC 1.14.14.1), CYP2U1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CYP2U1

Image



Immunofluorescent analysis of HeLa cells using CSB-PA768766LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.