



ISCA2 Antibody, FITC conjugated

Product Code	CSB-PA769773LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86U28
Immunogen	Recombinant Human Iron-sulfur cluster assembly 2 homolog, mitochondrial protein (9-154AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iron-sulfur cluster assembly 2 homolog, mitochondrial (HESB-like domain-containing protein 1), ISCA2, HBLD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ISCA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.