



RSPH3 Antibody, HRP conjugated

| | |
|---------------------|---|
| Product Code | CSB-PA769775LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q86UC2 |
| Immunogen | Recombinant Human Radial spoke head protein 3 homolog protein (1-275AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Radial spoke head protein 3 homolog (A-kinase anchor protein RSPH3) (Radial spoke head-like protein 2), RSPH3, RSHL2 RSP3 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | RSPH3 |