



RSPH3 Antibody, HRP conjugated

Product Code	CSB-PA769775LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UC2
Immunogen	Recombinant Human Radial spoke head protein 3 homolog protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Radial spoke head protein 3 homolog (A-kinase anchor protein RSPH3) (Radial spoke head-like protein 2), RSPH3, RSHL2 RSP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RSPH3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.