



RSPH3 Antibody, Biotin conjugated

Product Code	CSB-PA769775LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UC2
Immunogen	Recombinant Human Radial spoke head protein 3 homolog protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Radial spoke head protein 3 homolog (A-kinase anchor protein RSPH3) (Radial spoke head-like protein 2), RSPH3, RSHL2 RSP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RSPH3