




TBC1D26 Antibody

Product Code	CSB-PA769776LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UD7
Immunogen	Recombinant Human TBC1 domain family member 26 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 26, TBC1D26
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TBC1D26

Image



Western blot
 All lanes: TBC1D26 antibody at 3µg/ml + Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 29, 25, 10 kDa
 Observed band size: 29 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.