



LCA5 Antibody, Biotin conjugated

Product Code	CSB-PA769789LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VQ0
Immunogen	Recombinant Human Lebercilin protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lebercilin (Leber congenital amaurosis 5 protein), LCA5, C6orf152
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LCA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.