





OOSP2 Antibody

Product Code	CSB-PA769798LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WS3
Immunogen	Recombinant Human Oocyte-secreted protein 2 protein (74-158AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oocyte-secreted protein 2 (Placenta-specific 1-like protein) (Protein TMEM122), OOSP2, PLAC1L TMEM122
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OOSP2

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA769798LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA769798LA01HU at dilution of 1:100