



OOSP2 Antibody, Biotin conjugated

Product Code	CSB-PA769798LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WS3
Immunogen	Recombinant Human Oocyte-secreted protein 2 protein (74-158AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oocyte-secreted protein 2 (Placenta-specific 1-like protein) (Protein TMEM122), OOSP2, PLAC1L TMEM122
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OOSP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.