



NUDT13 Antibody, Biotin conjugated

Product Code	CSB-PA769800LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86X67
Immunogen	Recombinant Human Nucleoside diphosphate-linked moiety X motif 13 protein (32-263AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoside diphosphate-linked moiety X motif 13 (Nudix motif 13) (EC 3.-.-.-) (Protein KiSS-16), NUDT13
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NUDT13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.