



PPP1R3B Antibody, Biotin conjugated

Product Code	CSB-PA769802LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86XI6
Immunogen	Recombinant Human Protein phosphatase 1 regulatory subunit 3B protein (1-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1 regulatory subunit 3B (Hepatic glycogen-targeting protein phosphatase 1 regulatory subunit GL) (Protein phosphatase 1 regulatory subunit 4) (PP1 subunit R4) (Protein phosphatase 1 subunit GL) (PTG), PPP1R3B, PPP1R4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPP1R3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.