



CSMD2 Antibody, FITC conjugated

Product Code	CSB-PA770346LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z408
Immunogen	Recombinant Human CUB and sushi domain-containing protein 2 protein (1-225AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CUB and sushi domain-containing protein 2 (CUB and sushi multiple domains protein 2), CSMD2, KIAA1884
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CSMD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.