



FAM19A1 Antibody, Biotin conjugated

Product Code	CSB-PA770353LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z5A9
Immunogen	Recombinant Human Protein FAM19A1 protein (20-133AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM19A1 (Chemokine-like protein TAFA-1), FAM19A1, TAFA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM19A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.