




# HSD11B1L Antibody

<b>Product Code</b>	CSB-PA770354LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q7Z5J1
<b>Immunogen</b>	Recombinant Human Hydroxysteroid 11-beta-dehydrogenase 1-like protein (17-87AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Hydroxysteroid 11-beta-dehydrogenase 1-like protein (EC 1.1.1.-) (11-beta-hydroxysteroid dehydrogenase type 3) (11-DH3) (11-beta-HSD3) (Short chain dehydrogenase/reductase family 26C member 2) (Short-chain dehydrogenase/reductase 10), HSD11B1L, HSD3 SCDR10 SDR26C2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	HSD11B1L

## Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA770354LA01HU at dilution of 1:100

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.