



HSD11B1L Antibody, HRP conjugated

Product Code	CSB-PA770354LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z5J1
Immunogen	Recombinant Human Hydroxysteroid 11-beta-dehydrogenase 1-like protein (17-87AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydroxysteroid 11-beta-dehydrogenase 1-like protein (EC 1.1.1.-) (11-beta-hydroxysteroid dehydrogenase type 3) (11-DH3) (11-beta-HSD3) (Short chain dehydrogenase/reductase family 26C member 2) (Short-chain dehydrogenase/reductase 10), HSD11B1L, HSD3 SCDR10 SDR26C2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HSD11B1L