



HSD11B1L Antibody, HRP conjugated

Product Code	CSB-PA770354LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z5J1
Immunogen	Recombinant Human Hydroxysteroid 11-beta-dehydrogenase 1-like protein (17-87AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydroxysteroid 11-beta-dehydrogenase 1-like protein (EC 1.1.1.-) (11-beta-hydroxysteroid dehydrogenase type 3) (11-DH3) (11-beta-HSD3) (Short chain dehydrogenase/reductase family 26C member 2) (Short-chain dehydrogenase/reductase 10), HSD11B1L, HSD3 SCDR10 SDR26C2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HSD11B1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.