



SPATA19 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA770355LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q7Z5L4 |
| Immunogen | Recombinant Human Spermatogenesis-associated protein 19, mitochondrial protein (25-167AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Spermatogenesis-associated protein 19, mitochondrial (Spermatogenic cell-specific gene 1 protein) (Spergen-1), SPATA19, SPERGEN1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cancer |
| Target Names | SPATA19 |