



GALNT7 Antibody, Biotin conjugated

Product Code	CSB-PA773016LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86SF2
Immunogen	Recombinant Human N-acetylgalactosaminyltransferase 7 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acetylgalactosaminyltransferase 7 (EC 2.4.1.-) (Polypeptide GalNAc transferase 7) (GalNAc-T7) (pp-GaNTase 7) (Protein-UDP acetylgalactosaminyltransferase 7) (UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 7), GALNT7
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GALNT7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.