



# SMG6 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA773032LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q86US8
<b>Immunogen</b>	Recombinant Human Telomerase-binding protein EST1A protein (251-511AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Telomerase-binding protein EST1A (EC 3.1.-.-) (EST1-like protein A) (Ever shorter telomeres 1A) (Smg-6 homolog) (Telomerase subunit EST1A) (hSmg5/7a), SMG6, C17orf31 EST1A KIAA0732
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cancer
<b>Target Names</b>	SMG6
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.