



# LYZL2 Antibody, FITC conjugated

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA773594LC01HU                                                               |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| <b>Uniprot No.</b>         | Q7Z4W2                                                                           |
| <b>Immunogen</b>           | Recombinant Human Lysozyme-like protein 2 protein (63-194AA)                     |
| <b>Raised In</b>           | Rabbit                                                                           |
| <b>Species Reactivity</b>  | Human                                                                            |
| <b>Form</b>                | Liquid                                                                           |
| <b>Conjugate</b>           | FITC                                                                             |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| <b>Purification Method</b> | >95%, Protein G purified                                                         |
| <b>Isotype</b>             | IgG                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                       |
| <b>Alias</b>               | Lysozyme-like protein 2 (Lysozyme-2) (EC 3.2.1.17), LYZL2                        |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                             |
| <b>Research Area</b>       | Others                                                                           |
| <b>Target Names</b>        | LYZL2                                                                            |