




C5orf24 Antibody

Product Code	CSB-PA773601LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z6I8
Immunogen	Recombinant Human UPF0461 protein C5orf24 protein (7-188AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	UPF0461 protein C5orf24, C5orf24
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C5orf24

Image



Western blot

All lanes: C5orf24 antibody at 0.3µg/ml

Lane 1: Mouse gonadal tissue

Lane 2: K562 whole cell lysate

Lane 3: THP-1 whole cell lysate

Lane 4: U87 whole cell lysate

Lane 5: HepG2 whole cell lysate


Lane 6: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 21, 17 kDa

Observed band size: 17 kDa



Immunohistochemistry of paraffin-embedded

human kidney tissue using CSB-

PA773601LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.