





SLC46A3 Antibody

Product Code	CSB-PA773784LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z3Q1
Immunogen	Recombinant Human Solute carrier family 46 member 3 protein (26-73AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 46 member 3, SLC46A3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC46A3


Image




Western blot
All lanes: SLC46A3 antibody at 4µg/ml + Jurkat whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 52, 43 kDa
Observed band size: 52 kDa



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA773784LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA773784LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA773784LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.