



SLC46A3 Antibody, HRP conjugated

Product Code	CSB-PA773784LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z3Q1
Immunogen	Recombinant Human Solute carrier family 46 member 3 protein (26-73AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 46 member 3, SLC46A3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC46A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.