



GPR141 Antibody

Product Code	CSB-PA773798LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z602
Immunogen	Recombinant Human Probable G-protein coupled receptor 141 protein (153-183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 141 (G-protein coupled receptor PGR13), GPR141, PGR13
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR141
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.