





ARPIN Antibody

Product Code	CSB-PA773801LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z6K5
Immunogen	Recombinant Human Arpin protein (1-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arpin (Arp2/3 inhibition protein), ARPIN, C15orf38
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ARPIN


Image



Western blot
All lanes: ARPIN antibody at 0.5µg/ml + HepG2 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 25, 44 kDa
Observed band size: 25 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA773801LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA773801LA01HU at dilution of 1:100