



FAM131B Antibody

Product Code	CSB-PA774828LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86XD5
Immunogen	Recombinant Human Protein FAM131B protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM131B, FAM131B, KIAA0773
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM131B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.