




KRT73 Antibody

Product Code	CSB-PA774833LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86Y46
Immunogen	Recombinant Human Keratin, type II cytoskeletal 73 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type II cytoskeletal 73 (Cytokeratin-73) (CK-73) (Keratin-73) (K73) (Type II inner root sheath-specific keratin-K6irs3) (Type-II keratin Kb36), KRT73, K6IRS3 KB36 KRT6IRS3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KRT73

Image



Western blot
 All lanes: KRT73 antibody at 4µg/ml + HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 59, 42 kDa
 Observed band size: 59, 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.