



ERO1B Antibody, Biotin conjugated

Product Code	CSB-PA774834LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YB8
Immunogen	Recombinant Human ERO1-like protein beta protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ERO1-like protein beta (ERO1-L-beta) (EC 1.8.4.-) (Endoplasmic reticulum oxidoreductase beta) (Endoplasmic reticulum oxidoreductin-1-like protein B) (Oxidoreductin-1-L-beta), ERO1B, ERO1LB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ERO1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.