



PID1 Antibody, HRP conjugated

Product Code	CSB-PA800096LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z2X4
Immunogen	Recombinant Human PTB-containing, cubilin and LRP1-interacting protein (15-168AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PTB-containing, cubilin and LRP1-interacting protein (P-CLI1) (Phosphotyrosine interaction domain-containing protein 1) (Protein NYGGF4), PID1, NYGGF4 PCLI1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PID1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.