



PID1 Antibody, Biotin conjugated

Product Code	CSB-PA800096LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z2X4
Immunogen	Recombinant Human PTB-containing, cubilin and LRP1-interacting protein (15-168AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PTB-containing, cubilin and LRP1-interacting protein (P-CL11) (Phosphotyrosine interaction domain-containing protein 1) (Protein NYGGF4), PID1, NYGGF4 PCL11
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PID1