





RD3 Antibody

Product Code	CSB-PA800101LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z3Z2
Immunogen	Recombinant Human Protein RD3 protein (1-195AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein RD3 (Retinal degeneration protein 3), RD3, C1orf36
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RD3

Image



Western blot
All lanes: RD3 antibody at 0.5µg/ml
Lane 1: Mouse thymus tissue
Lane 2: Human placenta tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 23 kDa
Observed band size: 23 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA800101LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.