



RD3 Antibody, Biotin conjugated

Product Code	CSB-PA800101LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z3Z2
Immunogen	Recombinant Human Protein RD3 protein (1-195AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein RD3 (Retinal degeneration protein 3), RD3, C1orf36
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.