



SAMD7 Antibody, FITC conjugated

Product Code	CSB-PA801801LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z3H4
Immunogen	Recombinant Human Sterile alpha motif domain-containing protein 7 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sterile alpha motif domain-containing protein 7 (SAM domain-containing protein 7), SAMD7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SAMD7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.