



VN1R5 Antibody, Biotin conjugated

Product Code	CSB-PA801813LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z5H4
Immunogen	Recombinant Human Vomeronasal type-1 receptor 5 protein (185-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vomeronasal type-1 receptor 5 (G-protein coupled receptor GPCR26) (hGPCR26) (V1r-like receptor 5), VN1R5, V1RL5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	VN1R5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.