



RTN4RL1 Antibody, HRP conjugated

Product Code	CSB-PA803121LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UN2
Immunogen	Recombinant Human Reticulon-4 receptor-like 1 protein (201-420AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Reticulon-4 receptor-like 1 (Nogo receptor-like 2) (Nogo-66 receptor homolog 2) (Nogo-66 receptor-related protein 3) (NgR3), RTN4RL1, NGRH2 NGRL2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RTN4RL1