



RTN4RL1 Antibody, FITC conjugated

Product Code	CSB-PA803121LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UN2
Immunogen	Recombinant Human Reticulon-4 receptor-like 1 protein (201-420AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Reticulon-4 receptor-like 1 (Nogo receptor-like 2) (Nogo-66 receptor homolog 2) (Nogo-66 receptor-related protein 3) (NgR3), RTN4RL1, NGRH2 NGRL2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RTN4RL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.