



TTC30A Antibody, FITC conjugated

Product Code	CSB-PA803138LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WT1
Immunogen	Recombinant Human Tetratricopeptide repeat protein 30A protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 30A (TPR repeat protein 30A), TTC30A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TTC30A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.