



HORMAD1 Antibody, HRP conjugated

Product Code	CSB-PA803140LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86X24
Immunogen	Recombinant Human HORMA domain-containing protein 1 protein (201-394AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HORMA domain-containing protein 1 (Cancer/testis antigen 46) (CT46) (Newborn ovary HORMA protein), HORMAD1, NOHMA
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	HORMAD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.