



ZGRF1 Antibody, HRP conjugated

Product Code	CSB-PA803153LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YA3
Immunogen	Recombinant Human Protein ZGRF1 protein (3-302AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ZGRF1 (GRF-type zinc finger domain-containing protein 1), ZGRF1, C4orf21
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZGRF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.