



Txnip Antibody, Biotin conjugated

Product Code	CSB-PA803849LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8BG60
Immunogen	Recombinant Mouse Thioredoxin-interacting protein (1-397AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin-interacting protein (Vitamin D3 up-regulated protein 1), Txnip, Vdup1
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Txnip
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.