



Paox Antibody, Biotin conjugated

Product Code	CSB-PA804896LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8C0L6
Immunogen	Recombinant Mouse Peroxisomal N(1)-acetyl-spermine/spermidine oxidase protein (1-504AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal N(1)-acetyl-spermine/spermidine oxidase (EC 1.5.3.13) (Polyamine oxidase), Paox, Pao
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Paox
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.