



Paox Antibody, Biotin conjugated

Product Code	CSB-PA804896LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8C0L6
Immunogen	Recombinant Mouse Peroxisomal N(1)-acetyl-spermine/spermidine oxidase protein (1-504AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal N(1)-acetyl-spermine/spermidine oxidase (EC 1.5.3.13) (Polyamine oxidase), Paox, Pao
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Paox