



C1orf210 Antibody

Product Code	CSB-PA808993LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVY1
Immunogen	Recombinant Human Type III endosome membrane protein TEMP protein (51-113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Type III endosome membrane protein TEMP (TEMP), C1orf210
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C1orf210
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.