




UBR2 Antibody

Product Code	CSB-PA809001LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWW8
Immunogen	Recombinant Human E3 ubiquitin-protein ligase UBR2 protein (101-439AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase UBR2 (EC 2.3.2.27) (N-recognin-2) (RING-type E3 ubiquitin transferase UBR2) (Ubiquitin-protein ligase E3-alpha-2) (Ubiquitin-protein ligase E3-alpha-II), UBR2, C6orf133 KIAA0349
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBR2

Image



Western Blot

All lanes: UBR2 antibody at 3µg/ml + Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 201, 50, 67 kDa

Observed band size: 201 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.