



ZSCAN29 Antibody, Biotin conjugated

Product Code	CSB-PA809002LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWIY8
Immunogen	Recombinant Human Zinc finger and SCAN domain-containing protein 29 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger and SCAN domain-containing protein 29 (Zinc finger protein 690), ZSCAN29, ZNF690
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZSCAN29