





C11orf53 Antibody

Product Code	CSB-PA809008LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXP5
Immunogen	Recombinant Human Uncharacterized protein C11orf53 protein (13-236AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C11orf53, C11orf53
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C11orf53


Image



Western blot
 All lanes: C11orf53 antibody at 0.8µg/ml +
 Mouse gonadal tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 26 kDa
 Observed band size: 26 kDa



Immunohistochemistry of paraffin-embedded
 human colon cancer using CSB-
 PA809008LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA809008LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.